

Mucin-4

Catalog No: #AP85779



Package Size: #AP85779-1 50ug #AP85779-2 100ug #AP85779-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mucin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MUC-4,Ascites sialoglycoprotein,ASGP,Pancreatic adenocarcinoma mucin,Pre-sialomucin complex,pSMC,Sialomucin complex,Testis mucin,Smc
Accession No.	Uniprot ID:Q63661
Uniprot	Q63661
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Muc4

Note: This product is for in vitro research use only