

LETM1 domain-containing protein 1

Catalog No: #AP85781



Package Size: #AP85781-1 50ug #AP85781-2 100ug #AP85781-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LETM1 domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cervical cancer 1 proto-oncogene protein p40,Cervical cancer proto-oncogene 2 protein,HCCR-1,HCRR-2
Accession No.	Uniprot ID:Q6P1Q0
Uniprot	Q6P1Q0
GeneID	25875;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LETMD1

Note: This product is for in vitro research use only