

Neurotrophin-3

Catalog No: #AP85786



Package Size: #AP85786-1 50ug #AP85786-2 100ug #AP85786-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neurotrophin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NT-3,HDNF,Nerve growth factor 2,NGF-2,Neurotrophic factor,Ntf-3
Accession No.	Uniprot ID:P20181
Uniprot	P20181
GeneID	18205;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ntf3

Note: This product is for in vitro research use only