

# Myeloid differentiation primary response protein MyD88

Catalog No: #AP85790

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP85790-1 50ug #AP85790-2 100ug #AP85790-3 1mg

## Description

Product Name	Myeloid differentiation primary response protein MyD88
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	RCJMB04_14h3
Accession No.	Uniprot ID:A5HNF6
Uniprot	A5HNF6
GeneID	420420;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MYD88

Note: This product is for in vitro research use only