

Bile acid receptor

Catalog No: #AP85793



Package Size: #AP85793-1 50ug #AP85793-2 100ug #AP85793-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bile acid receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Farnesoid X-activated receptor,Farnesol receptor HRR-1,Nuclear receptor subfamily 1 group H member 4,Retinoid X receptor-interacting protein 14,RXR-interacting protein 14,Bar,Fxr,Rip14
Accession No.	Uniprot ID:Q60641
Uniprot	Q60641
GeneID	20186;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nr1h4

Note: This product is for in vitro research use only