

Nuclear factor interleukin-3-regulated protein

Catalog No: #AP85796



Package Size: #AP85796-1 50ug #AP85796-2 100ug #AP85796-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear factor interleukin-3-regulated protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	E4 promoter-binding protein 4,cE4BP4 protein,bZIP protein E4BP4,E4BP4
Accession No.	Uniprot ID:Q90Z72
Uniprot	Q90Z72
GeneID	395326;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NFIL3

Note: This product is for in vitro research use only