

Methylmalonyl-CoA mutase, mitochondrial

Catalog No: #AP85803



Package Size: #AP85803-1 50ug #AP85803-2 100ug #AP85803-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Methylmalonyl-CoA mutase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MCM,Methylmalonyl-CoA isomerase
Accession No.	Uniprot ID:P22033
Uniprot	P22033
GeneID	4594;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MUT

Note: This product is for in vitro research use only