

MOB kinase activator 1A

Catalog No: #AP85805



Package Size: #AP85805-1 50ug #AP85805-2 100ug #AP85805-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MOB kinase activator 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mob1 homolog 1B,Mps one binder kinase activator-like 1B,Mobk1b,Mobk11b
Accession No.	Uniprot ID:Q921Y0
Uniprot	Q921Y0
GeneID	232157;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mob1a

Note: This product is for in vitro research use only