

Dual specificity mitogen-activated protein kinase kinase 6

Catalog No: #AP85806

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85806-1 50ug #AP85806-2 100ug #AP85806-3 1mg

Description

Product Name	Dual specificity mitogen-activated protein kinase kinase 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAP kinase kinase 6,MAPK/ERK kinase 6,MEK 6,SAPKK3,Prkmk6,Sapkk3
Accession No.	Uniprot ID:P70236
Uniprot	P70236
GeneID	26399;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Map2k6

Note: This product is for in vitro research use only