

Canalicular multispecific organic anion transporter 2

Catalog No: #AP85809



Package Size: #AP85809-1 50ug #AP85809-2 100ug #AP85809-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Canalicular multispecific organic anion transporter 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ATP-binding cassette sub-family C member 3,MRP-like protein 2,MLP-2,Multidrug resistance-associated protein 3,Cmoat2,Mlp2,Mrp3
Accession No.	Uniprot ID:O88563
Uniprot	O88563
GeneID	140668;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abcc3

Note: This product is for in vitro research use only