

Neural cell adhesion molecule 2

Catalog No: #AP85810



Package Size: #AP85810-1 50ug #AP85810-2 100ug #AP85810-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neural cell adhesion molecule 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	N-CAM-2, Neural cell adhesion molecule RB-8, R4B12, Ocam, Rncam
Accession No.	Uniprot ID:O35136
Uniprot	O35136
GeneID	17968;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ncam2

Note: This product is for in vitro research use only