

# Myeloid differentiation primary response protein MyD88

Catalog No: #AP85816

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP85816-1 50ug #AP85816-2 100ug #AP85816-3 1mg

## Description

Product Name	Myeloid differentiation primary response protein MyD88
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Accession No.	Uniprot ID:A2TF48
Uniprot	A2TF48
GeneID	396646;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MYD88

Note: This product is for in vitro research use only