

Myogenin

Catalog No: #AP85817



Package Size: #AP85817-1 50ug #AP85817-2 100ug #AP85817-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myogenin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Class C basic helix-loop-helix protein 3,bHLHc3,Myogenic factor 4,Myf-4,BHLHC3,MYF4
Accession No.	Uniprot ID:P15173
Uniprot	P15173
GeneID	4656;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYOG

Note: This product is for in vitro research use only