

Metallothionein-2

Catalog No: #AP85825



Package Size: #AP85825-1 50ug #AP85825-2 100ug #AP85825-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Metallothionein-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MT-2,Metallothionein-II,MT-II,MT2A
Accession No.	Uniprot ID:P68301
Uniprot	P68301
GeneID	404070;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MT2

Note: This product is for in vitro research use only