

Neuromedin-S

Catalog No: #AP85827



Package Size: #AP85827-1 50ug #AP85827-2 100ug #AP85827-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neuromedin-S
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q5H8A3
Uniprot	Q5H8A3
GeneID	129521;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NMS

Note: This product is for in vitro research use only