

## Multidrug resistance-associated protein 9

Catalog No: #AP85830



Package Size: #AP85830-1 50ug #AP85830-2 100ug #AP85830-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Multidrug resistance-associated protein 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ATP-binding cassette sub-family C member 12,MRP9
Accession No.	Uniprot ID:Q96J65
Uniprot	Q96J65
GeneID	94160;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: ABCC12

Note: This product is for in vitro research use only