

Nucleoside diphosphate kinase A

Catalog No: #AP85832



Package Size: #AP85832-1 50ug #AP85832-2 100ug #AP85832-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nucleoside diphosphate kinase A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	NDK A, Metastasis inhibition factor NM23, Tumor metastatic process-associated protein
Accession No.	Uniprot ID:Q05982
Uniprot	Q05982
GeneID	191575;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nme1

Note: This product is for in vitro research use only