

## Protein NDRG1

Catalog No: #AP85833



Package Size: #AP85833-1 50ug #AP85833-2 100ug #AP85833-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein NDRG1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Differentiation-related gene 1 protein,DRG-1,N-myc downstream-regulated gene 1 protein,Nickel-specific induction protein Cap43,Reducing agents and tunicamycin-responsive protein,RTP,Rit42,CAP43,DRG1,RTP
Accession No.	Uniprot ID:Q92597
Uniprot	Q92597
GeneID	10397;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:NDRG1

Note: This product is for in vitro research use only