

Mitofusin-2

Catalog No: #AP85835

Package Size: #AP85835-1 50ug #AP85835-2 100ug #AP85835-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitofusin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Hypertension-related protein 1,Mitochondrial assembly regulatory factor,HSG protein,Transmembrane GTPase MFN2,Kiaa0214,Marf
Accession No.	Uniprot ID:Q80U63
Uniprot	Q80U63
GeneID	170731;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mfn2

Note: This product is for in vitro research use only