

Myogenic factor 6

Catalog No: #AP85838



Package Size: #AP85838-1 50ug #AP85838-2 100ug #AP85838-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myogenic factor 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Myf-6, Class C basic helix-loop-helix protein 4, bHLHc4, Muscle-specific regulatory factor 4, BHLHC4, MRF4
Accession No.	Uniprot ID:P23409
Uniprot	P23409
GeneID	4618;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MYF6

Note: This product is for in vitro research use only