

## Myosin-11

Catalog No: #AP85840



Package Size: #AP85840-1 50ug #AP85840-2 100ug #AP85840-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myosin-11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Myosin heavy chain 11, Myosin heavy chain, smooth muscle isoform, SMMHC, KIAA0866
Accession No.	Uniprot ID:P35749
Uniprot	P35749
GeneID	4629;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MYH11

Note: This product is for in vitro research use only