

## Myosin light chain 1/3, skeletal muscle isoform

Catalog No: #AP85841



Package Size: #AP85841-1 50ug #AP85841-2 100ug #AP85841-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myosin light chain 1/3, skeletal muscle isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MLC1/MLC3,Myosin light chain alkali 1/2,Myosin light chain A1/A2
Accession No.	Uniprot ID:A0JNJ5
Uniprot	A0JNJ5
GeneID	317657;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MYL1

Note: This product is for in vitro research use only