

Myosin regulatory light polypeptide 9

Catalog No: #AP85843



Package Size: #AP85843-1 50ug #AP85843-2 100ug #AP85843-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myosin regulatory light polypeptide 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Myosin regulatory light chain 2, smooth muscle isoform, Myosin regulatory light chain 9, MYRL2
Accession No.	Uniprot ID:Q5E9E2
Uniprot	Q5E9E2
GeneID	531505;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MYL9

Note: This product is for in vitro research use only