

Putative pro-MCH-like protein 1

Catalog No: #AP85844



Package Size: #AP85844-1 50ug #AP85844-2 100ug #AP85844-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Putative pro-MCH-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Pro-MCH variant,Pro-melanin-concentrating hormone-like protein 1
Accession No.	Uniprot ID:Q16048
Uniprot	Q16048
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PMCHL1

Note: This product is for in vitro research use only