

# Melanin-concentrating hormone receptor 1

Catalog No: #AP85845

Package Size: #AP85845-1 50ug #AP85845-2 100ug #AP85845-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Melanin-concentrating hormone receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	MCH receptor 1,G-protein coupled receptor 24,MCH-1R,MCH1R,GPR24
Accession No.	Uniprot ID:Q9MZ01
Uniprot	Q9MZ01
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MCHR1

Note: This product is for in vitro research use only