

Pituitary adenylate cyclase-activating polypeptide type I receptor

Catalog No: #AP85849

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85849-1 50ug #AP85849-2 100ug #AP85849-3 1mg

Description

Product Name	Pituitary adenylate cyclase-activating polypeptide type I receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PACAP type I receptor
Accession No.	Uniprot ID:P70205
Uniprot	P70205
GeneID	11517;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adcyap1r1

Note: This product is for in vitro research use only