

Oxysterols receptor LXR-beta

Catalog No: #AP85850



Package Size: #AP85850-1 50ug #AP85850-2 100ug #AP85850-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oxysterols receptor LXR-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Liver X receptor beta,Nuclear receptor subfamily 1 group H member 2,Orphan nuclear receptor OR-1,Ubiquitously-expressed nuclear receptor,UR,Lxrb
Accession No.	Uniprot ID:Q62755
Uniprot	Q62755
GeneID	58851;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nr1h2

Note: This product is for in vitro research use only