

Natural resistance-associated macrophage protein 2

Catalog No: #AP85853



Package Size: #AP85853-1 50ug #AP85853-2 100ug #AP85853-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Natural resistance-associated macrophage protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NRAMP 2,Divalent cation transporter 1,Divalent metal transporter 1,DMT-1,Solute carrier family 11 member 2,DCT1,DMT1,NRAMP2,OK/SW-cl.20
Accession No.	Uniprot ID:P49281
Uniprot	P49281
GeneID	4891;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC11A2

Note: This product is for in vitro research use only