

Arachidonate 5-lipoxygenase

Catalog No: #AP85855



Package Size: #AP85855-1 50ug #AP85855-2 100ug #AP85855-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Arachidonate 5-lipoxygenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	5-LO
Accession No.	Uniprot ID:P12527
Uniprot	P12527
GeneID	25290;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Alox5

Note: This product is for in vitro research use only