

Leucine-rich repeat serine/threonine-protein kinase 2

Catalog No: #AP85857



Package Size: #AP85857-1 50ug #AP85857-2 100ug #AP85857-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Leucine-rich repeat serine/threonine-protein kinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Dardarin,PARK8
Accession No.	Uniprot ID:Q5S007
Uniprot	Q5S007
GeneID	120892;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LRRK2

Note: This product is for in vitro research use only