

Neurogenic locus notch homolog protein 4

Catalog No: #AP85859



Package Size: #AP85859-1 50ug #AP85859-2 100ug #AP85859-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurogenic locus notch homolog protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Notch 4,Int-3,Int3
Accession No.	Uniprot ID:P31695
Uniprot	P31695
GeneID	18132;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Notch4

Note: This product is for in vitro research use only