

Mitogen-activated protein kinase 8

Catalog No: #AP85861



Package Size: #AP85861-1 50ug #AP85861-2 100ug #AP85861-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MAP kinase 8,SAPK gamma,Stress-activated protein kinase JNK1,c-Jun N-terminal kinase 1,p54 gamma,Jnk1,Prkm8
Accession No.	Uniprot ID:P49185
Uniprot	P49185
GeneID	116554;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mapk8

Note: This product is for in vitro research use only