

Matrix metalloproteinase-23

Catalog No: #AP85863



Package Size: #AP85863-1 50ug #AP85863-2 100ug #AP85863-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-23
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MMP-23, metalloprotease in the female reproductive tract, MIFR, Mifr
Accession No.	Uniprot ID:O88272
Uniprot	O88272
GeneID	94339;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mmp23

Note: This product is for in vitro research use only