

Lysyl oxidase homolog 4

Catalog No: #AP85864



Package Size: #AP85864-1 50ug #AP85864-2 100ug #AP85864-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lysyl oxidase homolog 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Lysyl oxidase-like protein 4
Accession No.	Uniprot ID:Q8MJ24
Uniprot	Q8MJ24
GeneID	281904;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LOXL4

Note: This product is for in vitro research use only