

Lipid phosphate phosphatase-related protein type 1

Catalog No: #AP85865



Package Size: #AP85865-1 50ug #AP85865-2 100ug #AP85865-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lipid phosphate phosphatase-related protein type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Plasticity-related gene 3 protein, PRG-3, Kiaa4247, Prg3
Accession No.	Uniprot ID:Q8BFZ2
Uniprot	Q8BFZ2
GeneID	272031;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Lppr1

Note: This product is for in vitro research use only