

Ly6/PLAUR domain-containing protein 3

Catalog No: #AP85867



Package Size: #AP85867-1 50ug #AP85867-2 100ug #AP85867-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ly6/PLAUR domain-containing protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GPI-anchored metastasis-associated protein C4.4A homolog,Matrigel-induced gene C4 protein,MIG-C4,C4.4A,UNQ491/PRO1007
Accession No.	Uniprot ID:O95274
Uniprot	O95274
GeneID	27076;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LYPD3

Note: This product is for in vitro research use only