

Nuclear receptor subfamily 2 group C member 2

Catalog No: #AP85868



Package Size: #AP85868-1 50ug #AP85868-2 100ug #AP85868-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor subfamily 2 group C member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Orphan nuclear receptor TAK1, Orphan nuclear receptor TR4, Testicular receptor 4, TAK1, TR4
Accession No.	Uniprot ID:P49116
Uniprot	P49116
GeneID	7182;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NR2C2

Note: This product is for in vitro research use only