

N-acetylglutamate synthase, mitochondrial

Catalog No: #AP85870



Package Size: #AP85870-1 50ug #AP85870-2 100ug #AP85870-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-acetylglutamate synthase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Amino-acid acetyltransferase
Accession No.	Uniprot ID:Q8N159
Uniprot	Q8N159
GeneID	162417;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NAGS

Note: This product is for in vitro research use only