

C-type mannose receptor 2

Catalog No: #AP85872

Package Size: #AP85872-1 50ug #AP85872-2 100ug #AP85872-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-type mannose receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C-type lectin domain family 13 member E, Endocytic receptor 180, Macrophage mannose receptor 2, Urokinase-type plasminogen activator receptor-associated protein, UPAR-associated protein, CD280, CLEC13E, ENDO180, KIAA0709, UPARAP
Accession No.	Uniprot ID: Q9UBG0
Uniprot	Q9UBG0
GeneID	9902;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MRC2

Note: This product is for in vitro research use only