

Neural cell adhesion molecule 1

Catalog No: #AP85874



Package Size: #AP85874-1 50ug #AP85874-2 100ug #AP85874-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neural cell adhesion molecule 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	N-CAM-1,CD56,Ncam
Accession No.	Uniprot ID:P13595
Uniprot	P13595
GeneID	17967;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ncam1

Note: This product is for in vitro research use only