

Ornithine aminotransferase, mitochondrial

Catalog No: #AP85876



Package Size: #AP85876-1 50ug #AP85876-2 100ug #AP85876-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ornithine aminotransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Ornithine--oxo-acid aminotransferase
Accession No.	Uniprot ID:Q3ZCF5
Uniprot	Q3ZCF5
GeneID	505323;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: OAT

Note: This product is for in vitro research use only