

Neurabin-2

Catalog No: #AP85877



Package Size: #AP85877-1 50ug #AP85877-2 100ug #AP85877-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurabin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Neurabin-II,Protein phosphatase 1 regulatory subunit 9B,Spinophilin
Accession No.	Uniprot ID:Q6R891
Uniprot	Q6R891
GeneID	217124;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ppp1r9b

Note: This product is for in vitro research use only