

Mannan-binding lectin serine protease 1

Catalog No: #AP85881



Package Size: #AP85881-1 50ug #AP85881-2 100ug #AP85881-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mannan-binding lectin serine protease 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Complement factor MASP-3, Complement-activating component of Ra-reactive factor, Mannose-binding lectin-associated serine protease 1, MASP-1, Mannose-binding protein-associated serine protease, Ra-reactive factor serine protease p100, RaRF, Serine protease 5, CRARF, CRARF1, PRSS5
Accession No.	Uniprot ID:P48740
Uniprot	P48740
GeneID	5648;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MASP1

Note: This product is for in vitro research use only