

Olfactory receptor 51E2

Catalog No: #AP85885



Package Size: #AP85885-1 50ug #AP85885-2 100ug #AP85885-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Olfactory receptor 51E2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HPRAJ,Olfactory receptor OR11-16,Prostate-specific G-protein coupled receptor,PSGR
Accession No.	Uniprot ID:Q9H255
Uniprot	Q9H255
GeneID	81285;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OR51E2

Note: This product is for in vitro research use only