

Lipoamide acyltransferase component of branched-chain alpha-keto acid dehydrogenase complex, mitochondrial

Catalog No: #AP85887

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85887-1 50ug #AP85887-2 100ug #AP85887-3 1mg

Description

Product Name	Lipoamide acyltransferase component of branched-chain alpha-keto acid dehydrogenase complex, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Branched-chain alpha-keto acid dehydrogenase complex component E2,BCKAD-E2,Dihydrolipoamide acetyltransferase component of branched-chain alpha-keto acid dehydrogenase complex,Dihydrolipoamide branched chain transacylase,Dihydrolipoalysine-residue (2-methylpropanoyl)transferase
Accession No.	Uniprot ID:P11181
Uniprot	P11181
GeneID	280759;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DBT

Note: This product is for in vitro research use only