

Nuclear receptor corepressor 2

Catalog No: #AP85889



Package Size: #AP85889-1 50ug #AP85889-2 100ug #AP85889-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor corepressor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	N-CoR2, CTG repeat protein 26, SMAP270, Silencing mediator of retinoic acid and thyroid hormone receptor, SMRT, T3 receptor-associating factor, TRAC, Thyroid-, retinoic-acid-receptor-associated corepressor, CTG26
Accession No.	Uniprot ID:Q9Y618
Uniprot	Q9Y618
GeneID	9612;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NCOR2

Note: This product is for in vitro research use only