

N-acylethanolamine-hydrolyzing acid amidase

Catalog No: #AP85893



Package Size: #AP85893-1 50ug #AP85893-2 100ug #AP85893-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-acylethanolamine-hydrolyzing acid amidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Acid ceramidase-like protein,N-acylsphingosine amidohydrolase-like,ASAH-like protein,ASAHL,PLT
Accession No.	Uniprot ID:Q02083
Uniprot	Q02083
GeneID	27163;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NAAA

Note: This product is for in vitro research use only