

Glutamate receptor ionotropic, NMDA 1

Catalog No: #AP85894



Package Size: #AP85894-1 50ug #AP85894-2 100ug #AP85894-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamate receptor ionotropic, NMDA 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GluN1, Glutamate [NMDA] receptor subunit zeta-1, N-methyl-D-aspartate receptor subunit NR1, NMD-R1, NMDAR1
Accession No.	Uniprot ID:Q05586
Uniprot	Q05586
GeneID	2902;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GRIN1

Note: This product is for in vitro research use only