

NAD(P)H-hydrate epimerase

Catalog No: #AP85895



Package Size: #AP85895-1 50ug #AP85895-2 100ug #AP85895-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAD(P)H-hydrate epimerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	AI-BP, Apolipoprotein A-I-interacting protein, AIBP
Accession No.	Uniprot ID:Q6QRN6
Uniprot	Q6QRN6
GeneID	404132;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NAXE

Note: This product is for in vitro research use only