

Motor neuron and pancreas homeobox protein 1

Catalog No: #AP85896



Package Size: #AP85896-1 50ug #AP85896-2 100ug #AP85896-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Motor neuron and pancreas homeobox protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Homeobox protein HB9,Hlxb9
Accession No.	Uniprot ID:Q9QZW9
Uniprot	Q9QZW9
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mnx1

Note: This product is for in vitro research use only